Docket No. 0546-1082 Appln. No. 10/550,945

APPENDIX:

The Appendix includes the following items:

- $oxed{igwedge}$ a Declaration
- □ Exhibits 1-6

IN THE U.S. PATENT AND TRADEMARK OFFICE

In re application of

Steve HEALD et al.

Conf. 4057

Application No. 10/550,945

Group 1651

Filed November 30, 2006

Examiner Irene Marx

PREPARATION OF VANILLIN FROM MICROBIAL TRANSFORMATION MEDIA BY EXTRACTION BY
MEANS SUPERCRITICAL FLUIDS OR GASES

DECLARATION

Assistant Commissioner for Patents

P.O. Box 1450

Alexandria, VA 22313-1450

Sir:

- I, Mylène DARRICAU, am a French citizen. I am employed by the French firm SAFISIS (a subsidiary of LESAFFRE), in the capacity of R&D Project Manager, and have held that position since January 5th, 2009. My background is a Master of Engineering in Biotechnology with 7 years of experience in Fermentation and Microbiology; I have a diploma from Engineering School of INSA in Toulouse (France), and I am In charge of Research and Development projects for SAFISIS in Biomolecules.
 - I, Evelyne FONCHY-PENOT, hereby declare that I am a French citizen. I am employed by the French firm LESAFFRE INTERNATIONAL (a subsidiary of LESAFFRE), in the capacity of R&D Project Manager, and have held that position since 2000. LESAFFRE INTERNATIONAL is a service company which supports the different subsidiaries of LESAFFRE particularly in Research projects.

My background is PhD in biotechnology with more than 10 years experience in fermentation and Microbiology; I have a diploma from Engineer School from ESII. In Marseille (France), and I am in charge of Research and Development projects for LESAFFRE INTERNATIONAL in Biomolecules and Bioprocess.

- We are familiar with the above-identified U.S. patent application, its prosecution before the United States Patent and Trademark Office, and the applied references of Rabenhorst et al. (U.S. Patent No. 6,133,003), Muheim et al. (U.S. Patent No. 6,235,507) and Makin (U.S. Patent No. 4,474.994).
- In order to demonstrate the patentability of the present invention, we are submitting the following observations.

We confirm hereby that we have conducted the following experiments.

We state that, since my company got reports from independent laboratories that compare DSM 9991 and DSM 9992 of H&R patent and conclude that these two strains cannot be distinguished from one another; experiments were conducted on the sole DSM 9992 strain, exactly as disclosed in the H&R patent Examples 1 and 3.

We, Mrs DARRICAU and Mrs FONCHY-PENOT, who are researchers at SAF-ISIS and LESAFFRE INTERNATIONAL respectively, have realized the works detailed in the enclosed report (Exhibit 1)

As a summary the conclusions are that:

Strains DSM 9992 from US Patent No. 6,133,003 and Zyl 926 from the SAF-ISIS present application do not react at all the same way when applying, the same process as disclosed in examples 1 and 3 of US Patent No. 6,133,003.

Moreover, the HPLC analysis controls show that the ferulic acid was accumulated in the wort so that no bioconversion occurred with any of these two strains.

This means that the process in the US Patent No. 6,133,003is inoperatively disclosed in order to get the alleged results regarding vanillin production.

As a consequence we state that the prior art was duly worked out with the intention to reduce to practice the working examples of US Patent No. 6,133,003, with the above mentioned absence of results.

We state that we made all reasonable due efforts, that a person of ordinary skill in the art should make, by taking all the information given in the specification such as for instance pH and oxygen maintenance, and also by using our general knowledge for example, to control the foaming that occurred during the fermentation by adding an antifoam.

Accordingly, we consider that, contrary to the DSM 9991 and 9992 strains of US Patent No. 6,133,003, the Zyl 926 microorganism, of the present application of SAF-ISIS, deposited within CABI as IMI 390106 for the purposes of patent procedures under the Budapest Treaty, provides unexpected results.

The latter are demonstrated notably by:

- -a consumption of ferulic acid with concomitant production of above 11 g/L of vanillin
- -a vanillin product obtained there from absent of odoriferous by-products and with none of guaïacol, vinyl guaïacol, eugenol and isoeugenol being present at more than 100ppm.

4. We declare that all statements made herein of our own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under \$1001 of Title 18 of the United States code and that such willful false statements may jeopardize the validity of the application or any patent issuing thereon.

Date June 29, 2010

Mylène DARRICAU R&D Product Manager SAF-ISIS Zone Artisanale Soustons 40140 – France

Jen?

Evelyne FONCHY-PENOT R&D Project Manager LESAFFRE INTERNATIONNAL Rue Gabriel Peri 147 Marcq en Baroeul 59700 – France



Report on the "Characterisation of Strains belonging to the Genus Amycolatopsis by 16S rRNA Gene Sequencing"

Professor M. Goodfellow

Background

Initially, the client wished to have the taxonomic relationships between the three Amycolatopsis strains determined by 16S rRNA gene sequencing. This study was designed to establish the inferred phylogenetic relationships between the cultures and their position (s) in the Amycolatopsis 16S rRNA gene tree.

Receipt of Strains

The package containing the strains arrived intact on May 20th (See photographs attached).

Materials and Methods

Cultivation of strains. The strains were grown in shake flasks (150 rpm) in ISP 2 broth for 48 hours at 37°C. Following incubation, biomass was harvested by centrifugation.

16S rRNA gene sequencing analyses. Chromosomal DNA was isolated and 16S rRNA PCR was done using well established procedure (Lane, 1991). PCR product was purified and direct sequencing of 16S rRNA genes of the three Amycolatopsis strains was carried out using BigDveTM version 3.1 Terminator Cycle Sequencing Kit (PE Applied Biosystems, MA). Samples were completed using 10 ng of DNA for every 100 bp of template DNA and 3.2 pmol of primer. Post-sequencing reactions were purified away from unincorporated dye terminators using Performa DTR spincolumns (EdgeBiosystems, MD). Sequences were determined using an ABI 3730 Sequencer. The resultant almost complete 16S rRNA gene sequences (1360 nucleotides) were aligned manually against corresponding sequences of the type strains of Amycolatopsis species, retrieved from GenBank, using PHYDIT software (Chun, 1995). A phylogenetic tree was inferred by using neighbour-joining (Saitou & Nei, 1987) tree-making algorithm from the PHYDIT program and an evolutionary distance matrix generated by using the distance model of Jukes & Cantor (1969), The robustness of the phyletic lines recovered in the Amycolatopsis tree were determined in a bootstrap analysis (Felsenstein, 1985) of the neighbour-joining dataset by using the SEQBOOT and CONSENSE options from the PHYLIP suite of programs (Felsenstein, 1993).

Results

It is apparent from Figures 1 to 3 that the 16S rRNA gene sequences obtained for the three Amycolatopsis strains are of high quality as there are no ambiguities in the conserved regions of the traces. It is clear from the 16S rRNA Amycolatopsis tree that the three strains belong to the Amycolatopsis methanolica subclade, the taxonomic integrity of the latter is underpinned by a 100% bootstrap value (Fig. 4). It is also apparent from the tree and from the associated nucleotide similarity and difference

matrix (Table 1) that the three strains share identical 16S rRNA gene sequences with one another.

Conclusions

- The three Amycolatopsis strains belong to the Amycolatopsis methanolica subclade in the 16S rRNA Amycolatopsis gene tree.
- The three Amycolatopsis strains shared identical 16S rRNA gene sequences with one another.

References

Chun, J. (1995). Computer-assisted classification and identification of actinomycetes. Ph.D. Thesis, University of Newcastle, Newcastle upon Tyne, UK.

Felsenstein, J. (1985). Confidence limits on phylogeny: An appropriate use of the bootstrap. *Evolution* 39, 7783-7791.

Felsenstein, J. (1993). PHYLIP- Phylogenetic Inference Package version 3.5.1. Department of Genetics, University of Washington, Seattle, USA.

Jukes, T.H. & Cantor, C.R. (1969). Evolution of protein molecules. In Mammalian Protein Metabolism, volume 3, pp. 21-132. Edited by H. N. Munro, Academic Press, New York.

Lane, D.J. (1991). 16S/238 rRNA Sequencing. In Nucleic Acid Techniques in Bacterial Systematics, pp. 115-148. Edited by E. Stackebrandt & M. Goodfellow. John Wiley & Sons, Chichester, New York.

Saitou, N. & Nei, M. (1987). The neighbour-joining method: a new method for reconstructing phylogenetic trees. *Molecular Biology and Evolution* 4, 406-425.

Attachments with the Report

Legends for Figures.

Figs. 1a-3b. Chromatograms of Samples Λ , B and C respectively (a- sequencing with forward primer; b- sequencing with reverse primer).

Fig. 4. Neighbour-joining tree based on nearly complete 16S rRNA gene sequences showing relationships between the three Amycolatopsis strains and between them and the type strains of Amycolatopsis species. Numbers at the nodes indicate the level of bootstrap support (%) based on neighbour-joining analysis of 1000 resampled datasets; only values above 50% are given. The scale bar indicates 0.02 substitutions per nucleotide position.

Photographs.

Photographs 1-5. Receipt of strains.

Tables.

Table 1. Nucleotide similarity and difference matrix with the three strains which share identical 16S rRNA gene sequences with one another.

WZ, Goodfellow
Professor Michael Goodfellow

CCFRA Technology Limited Station Road, Chipping Campden Gloucestershire GL55 6LD, United Kingdom

Tel: +44 (0)1386 842000 Fax: +44 (0)1386 842100 www.camoden.co.uk

Richard Tailon LeSaffre Safisis Zone Artisanale 40140 Soustons France





food and drink innovation

Report No: MB/REP/MMR/115732/00002-00004/2

Report Date: 30/06/09

CONFIDENTIAL CERTIFICATE OF ANALYSIS Microbiology Department

ncroniology Department MicroID

Notes

The results provided herein relate only to the items tested.

RiboPtint³, RiboGroup³ and RiboPrinterTM are trademarks of QualiconTM L.L.C., a subsidiary of B.L. DuPont de Nemours and Company, Wilmington, Delaware, USA.

Prepared/Checked By:

Authorised C. Mitchell
Research Technical Officer
Signatories:
Signatories:
Signatories:
C. Baylis
Methods Research Manager

Any opinions and interpretations contained in this report are outside the scope of UKAS accreditation.

This test report or certificate shall not be reproduced, except in full, without written approval of Campden BRL.

Confidentially notice. The information in this document is only intended for the individual or entity to whom it is addressed. It may contain privileged and confidential information that is extent from the indexed to the new of it you are not the intended recipitati, you must not copy, distribute or take any action to reliance on it. If you have received this document is every fewer needly copyed to lid intended by by telephone.

Report No: MB/REP/MMR/115732/00002-00004/2

Page 1 of 6

CCFRA Technology Lindted, Registered No. 1836922, Incorporated in England & Wales Registered Officer Station Road, Chipping Campden, Glos. GL55 6LD Part of Campden Bill Group

information envaluing from Campdes BRI is given after the oxercise of all reasonable care and skill in its considerion, preparation and issue, but its provided without liability in its application and use.

Summary

Three bacterial isolates were sent to Campden BRI from SAFISIS for strain comparison using the QualiconTM RiboPrinter® at the request of the client. The parcel was received on site unopened as shown in the images taken of the package in Appendix 1.

The isolates were stored and grown prior to analysis as detailed in the NDA supplied by the client. Standard techniques and reagents were used for the ribotyping. The enzymo PvuII was used for the analysis of the isolates as it is has been shown to differentiate between DuPont database entries for Streptomyces sp. The results of the work are noted in Tables below, with a pattern comparison shown in Figure I. A summary of the RiboPrinter methodology is outlined in Appendix 2.

CCFRA Code; Customer Code Sample Description: Comments:	MMR/115732/00002 Sample A Prozeu culture Stored at ambient temperature prior to analysis	Date of Receipt; Client presumptive identification	20/05/09 . Amycolatopsis sp
Analysis Type RiboPrinter®	Test Details Pruil	Identification None	Comments Isolate fell outside of 0.85 similarity needed for automatic identification against the DuPout identification database The isolate was placed in RiboGroup 114-3210-8-1.
RiboPrinter®	Strain comparison		Refer to Figure 1 for details


CCFRA Code: Customer Code Sample Description: Comments:	MMR/115732/00003 Sample B Freeze Dried Ampoule Stored at ambient temperature prior to analysis	Date of Receipt: Client presumptive identification	20/05/09 Amycolatopsis sp
Analysis Type RiboPrintea®	Test Details FyuII	Identification None	Comments Isolate fell outside of 0.85 similarity needed for automatic identification against the DuPont identification database. The isolate was placed in RiboGroup 114-3210-8-1.
RiboPrinter®	Strain comparison		Refer to Figure 1 for details


CCFRA Code: Customer Code Sample Description: Comments:	MMR/115732/00004 Sample C Freeze Dried Ampoule Stored at ambient temperature prior to analysis	Date of Receipt: Client presumptive identification	20/05/09 Amycolatopsis sp
Analysis Type RiboPrinter®	Test Details РуиН	Identification None	Comments Isolate fell outside of 0.85 similarity needed for automatic identification against the DuPont identification database. The isolate was placed in RiboGroup 114-3210-S-1.
RiboPrinter®	Strain comparison		Refer to Figure 1 for details

Strain comparison

A comparison of the patterns generated for the three isolates is shown below in Figure 1a and b.

Figure 1a: Strain comparison of isolates MB/115732/2-4 using MB/115732/2 as the reference isolate






All three isolates were placed in the same Ribogroup 114-3210-S-1. These results suggest that the isolates are indistinguishable from each other using PvuII.

These data readily demonstrate the ability of the RiboPrinter system to characterise and discriminate bacterial strains by analysis of RiboPrint patterns. This information will be available for reference on the archived database and as such can be used to authenticate, or compare to, fiese columns times time and time again.

These isolates will be kept for a standard period of 4 weeks after ribotyping after which they will be destroyed as outlined in the NDA supplied by the client.



The RiboPrinter system is an automated ribotyping process as follows:

Bacterial ceils are broken open to release the genomic DNA into solution. This DNA is digested ('cut up') using a restriction enzyme to generate fragments of different lengths, which are then separated out on the basis of size. Ribosomal DNA is detected using probes which give off light that is captured by a camera. The resulting image is converted into a RiboPrint pattern that can be characterised (finto RiboGroups) and where possible identified by proprietary software. RiboGroups are assigned to patterns in the following format - RiboGroup library - Machine member - batch number and sample position where the pattern occurred for the first time for example RIBO1 114-1002-S-1

As each sample is analysed by the software, it makes a decision about how closely the new sample pattern matches the existing patterns in the database. If the software finds a match, that is if it is mable to distinguish the new pattern from an existing pattern (i.e. >92% similarity between patterns), it assigns the new sample to that RiboGroup. All samples in a RiboGroup have the same RiboGroup code. Thus if the isolate has been previously ribotyped on this instrument, the instrument database will recognize the reproduced RiboFrint and allocate the RiboGroup code assigned in the first instance.

When an isolate is analysed on the machine for the first time, i.e. it does not match with any existing patterns, it will receive a new code i.e. the batch and sample number specific to that new batch run.

For automatic (DuPont) identification, the new sample information is compared to existing patterns within a supplied (DUP) database. If the software finds a match of >85% a DuPont identification will be reported.

If no automatic identification is given, the user can instruct the software to show the percentage similarities between existing patterns in the instrument sample datebase (Campden BRI database) and the new sample pattern. When >85% similarity is shown identification can be offered.



PHYDIT - Similarity analysis 1360 nucleotides analysed

Lower-left triangle contains [NT] Similarity.

Uppper-right triangle contains INTI Different/Total nucleotides.

Amycolatopsis fastidiosa DSM 43855T (AB184566)

The following table contains tab-delimited numbers. Copy and paste to MS Excel or Word.

Amycolatopsis Sample A Amycolatopsis Sample A Amycolatopsis Sample B 100 Amycolatopsis Sample C 100 Amycolatopsis methanolica DSM 44096T (AJ249135) 98.88 Amycolatopsis thermoflava DSM 44574T (AF052390) 98 88 Amycolatopsis jejuensis NRRL B-24427T (DQ000200) 94.25 Amycolatopsis sulphurea DSM 46092T (AJ293756) 95.37 Amycolatopsis alba DSM 44262T (AF051340) 94 76 Amycolatopsis coloradensis DSM 44225T (AJ293753) 94.47 Amycolatopsis regifaucium DSM 45072T (AY129760) 94.84 Amycolatopsis azurea DSM 43854T (AJ400709) 94.69 94.69 Amycolatopsis orientalis DSM 40040T (AJ400711) Amycolatopsis keratiniphila subsp. keratiniphila DSM 4440T (AJ278496) 94.69 94.69 Amycolatopsis keratiniphila subsp. nogabecina DSM 445T (AJ508238) Amycolatopsis japonica DSM 44213T (AJ508236) 94 54 Amycolatopsis decaplanina DSM 44594T (AJ508237) 94.77 Amycolatopsis minnesotensis NRRL B-24435 (DQ076482) 94 84 Amycolatopsis nigrescens DSM 44992T (DQ486888) 95.44 Amycolatopsis lexingtonensis DSM 44653T (AY183358) 94 61 Amycolatopsis pretoriensis DSM 44654T (AY183356) 94.54 94 69 Amycolatopsis rifamycinica DSM 46095T (AY083603) Amycolatopsis kentuckyensis DSM 44652T (AY183357) 94.54 Amycolatopsis balhimycina DSM 44591T (AJ508239) 94.47 Amycolatopsis vancoresmycina DSM 44592T (AJ508240) 94.84 94 61 Amycolatopsis plumensis DSM 44776T (AY262825) 94 61 Amycolatopsis tolypomycina DSM 44544T (AJ293757) Amycolatopsis mediterranei DSM 43304T (AJ293754) 94.84 Amycolatopsis australiensis DSM 44671T (AY129753) 95 14 Amycolatopsis saalfeldensis DSM 44993T (DQ792500) 95.07 Amycolatopsis echigonensis JCM 21831T (AB248535) 94.84 Amycolatopsis niigatensis JCM 21832T (AB248537) 95.06 Amycolatopsis benzoatilytica DSM 43387T (AY957506) 93.39 94.9 Amycolatopsis albidoflavus DSM 44639T (AJ252832) Amycolatopsis halotolerans NRRL B-24428T (DQ000196) 95.21 Amycolatopsis rubida DSM44637T (AF222022) 94 51 Amycolatopsis sacchari DSM 44468T (AF223354) 95.29 Amycolatopsis taiwanensis NBRC 102103T (DQ160215) 94 99 Amycolatopsis marina NBRC 104263T (EU329845) 95.96 Amycolatopsis palatopharyngis JCM 12460T (AF479268) 96.19

Amycolatopsis Sample B 0/1338	Amycolatopsis Sample C 0/1338 0/1338
100	
98.88	98.88
98.88	98.88
94.25	94.25
95.37	95.37
94.76	94.76
94,47	94.47
94.84	94.84
94.69	94.69
94.69	94.69
94.69	94.69
94.69	94.69
94.54	94.54
94.77	94.77
94.84	94.84
95.44	95.44
94.61	94.61
94.54	94.54
94.69	94.69
94.54	94.54
94.47	94.47
94.84	94.84
94.61	94.61
94.61	94.61
94.84	94.84
95.14	95.14
95.07	95.07
94.84	94.84
95.06	95.06
93.39	93.39
94.9	94.9
95.21	95.21
94.51	94.51
95.29	95.29
94.99	94.99
95.96	95.96
96.19	96.19

96.19 92.38

Amycolatopsis methanolica DSM 44096T (AJ249135)

15/1338

15/1338 15/1338

99.85 94.32

94.99

94.99 94.69

94.47

94.39

94.62 94.54

94.77

94.61

94.47 94.92

94.54

95.74

94.61

94.54

94.62

94.47

94.39

94.47

94.17 94.24

94.76

94.99

95.22

94.54

94.76

92.93

94.75

94.91

94.91

95.74

95,44

96.34

96.49 92.38 Amycolatopsis thermoflava DSM 44574T (AF052390)

15/1337

15/1337

15/1337 2/1337

-

94.24

94.99

94.61

94.39 94.31

94.54

94.47

94.69

94.54 94.39

94.84

94.47

95.66 94.54

94.47

94.54

94.39

94.32

94.39

94.09

94.16 94.68

94,91

95.14

94.46

94.68

92.85

94.67

94.83

94.73

95.66

95.36

96.26 96.41

Amycolatopsis jejuensis NRRL B-24427T (DQ000200)

77/1338

77/1338 77/1338

76/1338

77/1337

97.83 96.86

96.64 97.08

97.01

96.79 96.64

96.86

97.01 96.79

96.71

96.71

96.56

96.64 96.64

96.41

96.19

96.64

96.34

96.78

96.48

97.01

97.16

97.23

97.16

95.29

97

96.86 96.77

96.64

93.88

95.44

95.52

Amycolatopsis sulphurea DSM 46092T (AJ293756)	Amycolatopsis alba DSM 44262T (AF051340)
62/1338	70/1337
62/1338	70/1337
62/1338	70/1337
67/1338	71/1337
67/1337	72/1336
29/1339	42/1337
414	45/1337
96.63	
96.49	
97.38	
97.01	
96.94	
96.79	
97.01	
96.94	
97.31	
97.01	
97.09	
96.63	
96.56	
96.56	
96.34	
96.26	
96.71	
96.11	
96.41	
96.41	
97.09	
97.16	
97.31	
97.68	
95.67	
97.68	
97.68	
96.92	
97.23	
93.58	
95.52	
95.67	
91.04	91.1

Amycolatopsis coloradensis DSM 44225T (AJ293753)

74/1338

74/1338

74/1338

74/1338

75/1337 45/1338

47/1338

10/1337

98.5

98.58

98.8 98.65

98.65

98.65

98.58 97.61

97.38

97.83

97.91

97.83

97.61

97.68 97.46

97.08

97.53

97.68

97,46

97.91

97.08

97.23

95.59

97.45

97.23

97,14

95.89

93.87

95.07

95.37 90.81

```
Amycolatopsis regifaucium DSM 45072T (AY129760)
69/1337
69/1337
75/1337
```

76/1336 39/1337 35/1337

23/1336 20/1337

98.88 99.1

98.88 99.25

99.48 99.1

97,46

97.01

97.46 97.38

97.38

97.16

97.01

97.46

97.01

97.38 97.3

97.83

97.23

97.83

98.13

96.43

98.05

98.13

98.13

97.29

95.96

93.34

95.29 95.29

Amycolatopsis azurea DSM 43854T (AJ400709)	Amycolatopsis orientalis DSM 40040T (AJ4007	11)
71/1338	71/1338	
71/1338	71/1338	
71/1338	71/1338	
72/1338	73/1338	
73/1337	74/1337	
40/1338	43/1338	
40/1338	41/1338	
18/1337	19/1337	
19/1338	16/1338	
15/1337	12/1337	
	6/1338	
99.58		20.00
99.18		99.33 99.63
99.48		
99.		99.18
99.18		99.33
97.16 97.3		97.31 97.31
97.3		97.83
98.89		97.83 97.76
98.00		97.76 97.61
97.68		97.38
97.68		97.38
97.53		97.38
96.93		96.86
97.46		97.38
97.75		97.46
97.76		97.76
97.76		97.46
97.6		97.38
97.9		97.68
96.43		96.35
97.9		97.83
97.75		97.68
97.67		97.44
96.26		96.26
93.87		93.65
95.74		95.59
96.04		95.81
91.63		91.63

```
Amycolatopsis keratiniphila subsp. keratiniphila DSM 4440T (AJ278496)
```

71/1338

71/1338

71/1338

70/1338

71/1337 45/1338

43/1338

21/1337

18/1338

15/1337 11/1338

9/1338

99.7

99.25 99.33

97.31

97.53

97.83

97.76

97.46 97.23

97.23

97.46

96.86

97.31

97.38

97.61

97.38 97.76

97.9

96.05 97.98

97.75

97,89

96.49

93.87

95.89 96.11

Amycolatopsis keratiniphila subsp. nogabecina DSM 445T (AJ508238) 71/1337 71/1337 72/1337 73/1336 42/1337 40/1337

99.33 99.48 97.31 97.38 97.75 97.68 97.53 97.31 97.31 97.46 96.78 97.31 97.38 97.76 97.38 97.83 98.13 96.43 98.05 97.98 97.89 96.41 93.79 95.74 95.96 91.77

18/1337 10/1336 7/1337 5/1337

4/1337

Amycolatopsis japonica DSM 44213T (AJ50823	36)	Amycolatopsis decaplanina DSM 44594T (AJ508	3237)
73/1338		70/1338	
73/1338		70/1338	
73/1338		70/1338	
74/1338		68/1338	
75/1337		69/1337	
40/1338		43/1338	
41/1338		36/1338	
23/1337		22/1337	
18/1338		19/1338	
7/1337		12/1337	
12/1338		11/1338	
11/1338		9/1338	
10/1338		9/1338	
9/1337		7/1337	
		11/1338	
	9.18		07.00
	7.23		97.09
	7.16		97.23
	7.83		97.38 97.31
	7.61		97.31
	7.38		97.08
	7.23		96.94
	7.31		97.09
	6.93		96.48
	7.46		96.93
	7.46		97.31
	7.68		97.98
	7.38		97.09
	7.53		97.68
	7.75		98.13
	95.9		96.12
9	7.75		97.98
	97.6		97.98
9	7.37		97.67
9	5.96		96.19
	3.57		93.57
9	5.52		95.89
	5.59		95.81
9	1.11		91.41

Amycolatopsis minnesotensis NRRL B-24435 (DQ076482)

69/1338

69/1338

69/1338

73/1338

74/1337 44/1339

40/1339

35/1337

32/1338 34/1337

38/1338

36/1338 36/1338

36/1337 37/1338

39/1338

98.43 96.71

96.79 96.64

96.41 96.56

97.38

96.26

96.56 96.63

96.94

97.76 96.41

96.56

95.14 96.78

96.56

96.24 96.86

93.13

95.96 96.11

Amycolatopsis nigrescens DSM 44992T (DQ486888)

61/1338

61/1338

61/1338

57/1338 58/1337

44/1339

33/1337 35/1338

40/1337

36/1338 36/1338 33/1338

35/1337 38/1338

37/1338 21/1339

21/11

96.71

96.79 96.79

96.56

96.49 96.64

96.04

96.41

96.78 97.09

97.09

96.49

96.64 94.84

96.78

96.64 96.7

97.53 93.5

96.64

96.64 91.78

```
Amycolatopsis lexingtonensis DSM 44653T (AY183358) 72/1337 72/1337
```

72/1337 72/1337 73/1336

46/1337 45/1337

29/1336 29/1337

34/1336 25/1337

29/1337 29/1337 30/1336

29/1337 35/1337

44/1337 44/1337

99.93

99.63 99.4

99.25 98.95

98.8 99.1

99.1

98.28 98.5 97.6

97.75 96.65 97.97

97.75

97.59

96.34 94.32

95.36

95.51 91.47

```
Amycolatopsis pretoriensis DSM 44654T (AY183356) 73/1338 73/1338
```

73/1338 73/1338 73/1338 74/1337

45/1339 46/1339 28/1337

28/1338 35/1337 26/1338

30/1338 30/1338 31/1337

30/1338 36/1338 43/1339

43/1339 1/1337

__

99.7

99.48 99.33

99.03

98.88 99.18

99.18 98.36

98.58

97.53 97.68

96.58 97.9

97.68

97.52 96.41

94.4

95.29 95.44 91.41 Amycolatopsis rifamycinica DSM 46095T (AY083603)

71/1338

71/1338

71/1338

72/1338 73/1337

45/1338

46/1338

27/1337 29/1338

35/1337 28/1338

32/1338

34/1338

33/1337 32/1338

36/1338

45/1338

43/1338 5/1337

4/1338

99.63 99.18

99.03

98.65 98.95

99.33

98.36

98.43

97.68

97.83

96.43 97.9

97.83

97.52

96.34 94.32

95.37

95.37

```
Amycolatopsis kentuckyensis DSM 44652T (AY183357)
73/1337
73/1337
73/1337
74/1337
75/1336
48/1337
49/1337
30/1336
32/1337
38/1336
31/1337
35/1337
37/1337
36/1336
35/1337
39/1337
48/1337
46/1337
```

8/1336 7/1337 5/1337

> 98.95 98.58 98.88 99.25 98.28 98.2 97,46 97.6 96.2 97.67 97.75 97,44 96.04 94.02 95.21 95.21

91.1

```
Amycolatopsis balhimycina DSM 44591T (AJ508239)
74/1338
74/1338
74/1338
```

75/1338 76/1337 51/1338

50/1338

31/1338 40/1337 31/1338

35/1338 37/1338 36/1337

37/1338 41/1338 46/1338 47/1338

10/1337 9/1338

11/1338 12/1337

~~~

98.65 98.35

98.65 99.18 98.28

98.36 97.16

97.31 96.2 97.53

97.46 97.29

95.96 94.17

94.92 95.07

Amycolatopsis vancoresmycina DSM 44592T (AJ508240)

69/1338

69/1338

69/1338

74/1338

75/1337

45/1338

44/1338

31/1337 34/1338

34/1337

33/1338

35/1338

34/1338

36/1338

39/1338

35/1338

45/1338 14/1337

13/1338

13/1338

14/1337 18/1338

---

98.65

98.95

98.65

98.36

98.13 97.91

98.05

96.66

98.13

98.2

97.59

96.26

94.17

95.29

95.29

```
Amycolatopsis plumensis DSM 44776T (AY262825) 72/1337 72/1337
```

72/1337

78/1337 79/1336

49/1337 52/1337

38/1336 39/1337

40/1336

42/1337 42/1337

43/1336 41/1337 47/1337

50/1337 53/1337 16/1336

15/1337 18/1337 19/1336

22/1337 18/1337

---

99.1

98.28 97.83 97.46

97.38

97.45 96.27

97.45

97.6

96.84 95.74

93.79

94.76 94.99

```
Amycolatopsis tolypomycina DSM 44544T (AJ293757)
72/1337
72/1337
72/1337
77/1337
78/1336
43/1337
48/1337
32/1336
33/1337
35/1336
34/1337
35/1337
36/1337
36/1336
34/1337
41/1337
46/1337
```

48/1337 12/1336 11/1337 14/1337 15/1336 18/1337

> 14/1337 12/1336

\*\*\*

98.58 98.2 97.91 97.6 97.75 96.58 97.82 97.9 97.14 95.89 94.02 95.06 95.36

```
Amycolatopsis mediterranei DSM 43304T (AJ293754)
69/1336
69/1336
69/1336
70/1336
71/1335
47/1336
48/1336
28/1335
31/1336
36/1335
30/1336
34/1336
35/1336
35/1335
34/1336
36/1336
45/1336
43/1336
12/1335
11/1336
9/1336
```

23/1335 19/1335

10/1335 11/1336 18/1336

> 98.5 98.35 97.38 97.53 96.12 97.6 97.68 97.36 96.11 94.01 95.36 95.36 91.39

```
Amycolatopsis australiensis DSM 44671T (AY129753) 65/1338
```

65/1338

65/1338 67/1338

68/1337 40/1338

39/1338 34/1337

34/1338 29/1337

30/1338 30/1338 32/1338

30/1337 31/1338

27/1338 41/1338

39/1338 23/1337

22/1338

23/1337 23/1338 22/1338

29/1337 24/1337 20/1336

\_\_\_

97.91 97.46

97.75

96.35 97.98

97.9

97.52 96.71

94.25

95.96 95.81

91.41

Amycolatopsis saalfeldensis DSM 44993T (DQ792500)

66/1338

66/1338

66/1338

64/1338 65/1337

38/1338 38/1338

29/1337 28/1338

37/1337 30/1338

34/1338 35/1338

35/1337

35/1338 39/1338

30/1338 27/1338

20/1337

19/1338 21/1338

24/1337

22/1338

25/1338

34/1337 28/1337

22/1336 28/1338

97.16

97.31 95.74

97.53 97.31

97.22

97.01 94.39

95.74

95.89 91.48

```
Amycolatopsis echigonensis JCM 21831T (AB248535)
```

69/1337

69/1337

69/1337

73/1337 74/1336

37/1338

36/1338

35/1336 39/1337

29/1336 32/1337

35/1337

30/1337 29/1336

33/1337 31/1337

48/1338 47/1339

32/1336 33/1338

31/1337

34/1336 38/1337 28/1337

35/1336

32/1336 35/1335

34/1337 38/1337

99.4 97.19

98.95 99.1 98.5

96.19 94.02

95.44

95.29 90.88

```
Amycolatopsis niigatensis JCM 21832T (AB248537)
66/1336
66/1336
66/1336
70/1336
71/1335
38/1337
31/1337
33/1335
37/1336
25/1335
28/1336
31/1336
28/1336
25/1335
30/1336
25/1336
46/1337
45/1338
30/1335
31/1337
29/1336
32/1335
36/1336
26/1336
34/1335
```

30/1335 33/1334 30/1336 36/1336 8/1338

> 97.57 99.18 99.55 98.42 96.18 93.94 95.43 95.21

91.25

```
Amycolatopsis benzoatilytica DSM 43387T (AY957506)
87/1316
87/1316
87/1316
93/1316
94/1315
62/1316
57/1316
55/1315
58/1316
47/1315
47/1316
48/1316
52/1316
47/1315
54/1316
51/1316
64/1316
68/1317
44/1315
45/1316
47/1316
50/1315
```

50/1316 44/1316 49/1315 45/1315 51/1314 48/1316 56/1316 37/1316 32/1316

> 97.41 97.42 96.64 94.53 92.17 93.54 93.62

90.05

```
Amycolatopsis albidoflavus DSM 44639T (AJ252832)
68/1334
68/1334
68/1334
70/1334
71/1333
40/1335
31/1335
32/1333
34/1334
26/1333
28/1334
29/1334
27/1334
26/1333
30/1334
```

27/1334 43/1335 43/1336 27/1333 28/1335 28/1334 31/1333 33/1334 25/1334 34/1333 29/1333 32/1332 27/1334 33/1334

14/1336 11/1336 34/1314

> 99.18 98.5 96.33 93.93 95.35 95.28 91.32

```
Amycolatopsis halotolerans NRRL B-24428T (DQ000196)
64/1336
64/1336
64/1336
68/1336
69/1335
42/1337
31/1337
33/1335
37/1336
25/1335
30/1336
31/1336
30/1336
27/1335
32/1336
27/1336
46/1337
45/1338
30/1335
31/1337
29/1336
30/1335
34/1336
24/1336
32/1335
28/1335
31/1334
28/1336
36/1336
12/1338
```

6/1338 34/1316 11/1336

> 98.42 96.18 93.72 95.28 95.21 91.25

| 73/1330 | 63/1338 |       |
|---------|---------|-------|
| 73/1330 | 63/1338 |       |
| 73/1330 | 63/1338 |       |
| 69/1330 | 57/1338 |       |
| 70/1329 | 58/1337 |       |
| 43/1331 | 45/1338 |       |
| 41/1331 | 37/1338 |       |
| 37/1329 | 51/1337 |       |
| 38/1330 | 55/1338 |       |
| 36/1329 | 54/1337 |       |
| 31/1330 | 50/1338 |       |
| 34/1330 | 50/1338 |       |
| 28/1330 | 47/1338 |       |
| 28/1329 | 48/1337 |       |
| 35/1330 | 54/1338 |       |
| 31/1330 | 51/1338 |       |
| 50/1331 | 42/1338 |       |
| 44/1332 | 33/1338 |       |
| 32/1329 | 49/1337 |       |
| 33/1331 | 48/1338 |       |
| 33/1330 | 49/1338 |       |
| 34/1330 | 53/1337 |       |
| 36/1330 | 54/1338 |       |
| 32/1330 | 50/1338 |       |
| 42/1329 | 57/1337 |       |
| 38/1329 | 55/1337 |       |
| 35/1328 | 52/1336 |       |
| 33/1330 | 44/1338 |       |
| 37/1330 | 40/1338 |       |
| 20/1332 | 51/1337 |       |
| 21/1332 | 51/1336 |       |
| 44/1310 | 72/1316 |       |
| 20/1330 | 49/1334 |       |
| 21/1332 | 51/1336 |       |
| ~~~     | 46/1330 |       |
|         | 96.54   |       |
|         | 93.99   | 94.77 |
|         | 95.79   | 96.26 |
|         | 95.64   | 96.41 |
|         | 90.98   | 91.85 |

Amycolatopsis rubida DSM44637T (AF222022) Amycolatopsis sacchari DSM 44468T (AF223354)

```
Amycolatopsis taiwanensis NBRC 102103T (DQ160215)
67/1338
67/1338
67/1338
61/1338
62/1337
82/1339
86/1339
79/1337
82/1338
89/1337
82/1338
85/1338
82/1338
83/1337
86/1338
86/1338
92/1339
87/1339
76/1337
75/1339
76/1338
80/1337
78/1338
78/1338
83/1337
80/1337
80/1336
77/1338
75/1338
80/1338
81/1337
103/1316
81/1335
84/1337
```

80/1331 70/1338

> 94.1 94.02

> 89.92

Amycolatopsis marina NBRC 104263T (EU329845) 54/1338 54/1338 49/1338

50/1337 61/1338

54/1338

60/1338 65/1337

66/1338 63/1337 57/1338 59/1338

55/1338 57/1337 60/1338 55/1338

54/1338 45/1338 62/1337 63/1338 62/1338 64/1337

68/1338 63/1338 70/1337 66/1337 62/1336 54/1338 57/1338

61/1337 61/1336 85/1316 62/1334

63/1336 56/1330 50/1338

79/1338

99.4 91.55 Amycolatopsis palatopharyngis JCM 12460T (AF479268)

51/1338

51/1338

51/1338

47/1338 48/1337

60/1338 58/1338

61/1337

62/1338

63/1337

53/1338 56/1338

52/1338 54/1337 59/1338

59/1338 56/1338 52/1338

45/1338 60/1337 61/1338

62/1338 64/1337 66/1338

63/1338 67/1337 62/1337 62/1336

56/1338 55/1338 63/1337

63/1337 64/1336

84/1316 63/1334 64/1336

58/1330 48/1338

80/1338 8/1342 Amycolatopsis fastidiosa DSM 43855T (AB184566)

102/1338

102/1338

102/1338

102/1338 103/1337

128/1339

120/1339

123/1338

112/1338 112/1338

113/1338

119/1338 115/1338

120/1339 110/1339 114/1337

115/1339 115/1338 119/1337

119/1337 119/1338 116/1338

123/1337 119/1337 115/1336

115/1338 114/1338 122/1338

117/1337 131/1316 116/1336

116/1336 117/1337 120/1331 109/1338

135/1339 113/1338

110/1338

---

# Exhibit 4

#### >a1-27F

#### >a1-1492R

## >41-1492R Inv compl

#### Sea A

COCTA GOTTETIGGE GACTECCA COTTETICATION COTTACTA ACCA THANGGE COCCCCCCCCCGGGGGTACGG
COCCA AGCTA AACTA AAGCAA TEACACGCCCCCCCCACAGODY SACCA TO TIGGATHA ATTRATICA AC
OPGANAMA CITACCTOPS TEACA TEACACCCACAGODY SACCA TO TEACAGCC
OPGANAMA CITACCTOPS TEACACTOP THAT ACCACACAGODY SACCATOR TO THAT ACCACACAGODY ATTRATICA CACACAGODY ATTRATICA CACACAGODY ATTRATICA ACCACAGODY ATTRATICA CACACAGODY ATTRATICACAGODY ATTRATICACAGODY ATTRATICACAGODY ATTRATICACAGODY A

## >b2-27F

## >b2-1492R

#### >b2-1492R

A HOMPORGICT TRACTIC OBBANT TO LAGTICGATE ADDITIONATION TO GETTAGE TO OBTANGOGGAGAC I CGAPATENT TOUT OF A CO GETA ANTONIO AND THE PROBABILITY AND ADDITIONATION OF GETTAGE AND THE PROBABILITY AND ADDITIONATION OF A PROGAGGABOTH AND ADDITIONATION OF ADDITIONATIO


## Seq B

#### 561-279

## >c1-1492R

## >cl-1492P

## Seq C







## **EXTERNAL MEMORANDUM**

Re : Report of production trials according to US 6.133,003 to Rabenhorst et al.

## Table of contents

| 1. Goal                                                                                                                                                       |     |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------|-----|
| 2. Material and Methods                                                                                                                                       | 2   |
| 2.1 Cultures work out.       2         2.1.1 Inoculation chain.       2         2.1.2 Cultures in fermenters (max. working volume 15L).       3               |     |
| 2.2 Cultures follow up.       5         2.2.1 Culture parameters.       5         2.2.2 Bactenological follow up.       5         2.2.3 HPLC Analyses       6 |     |
| 3. Results                                                                                                                                                    |     |
| 3.1 Culture follow up.       6         3.1 1 DSM 9992.       6         3.1 2 Zyl 926 (Asp L-li 2).       8         3.1.3 Cultures comparison       9          |     |
| 3.2 Analytical Results                                                                                                                                        |     |
| A Complicing                                                                                                                                                  | 3 . |



#### 1. Goal

The goal is to compare two strains:

- Arnycolatopsis sp. DSM 9992
- Amycolatopsis sp. Zyl 926 (referred to also as Asp L-ii 2).

when cultivated in identical conditions, namely identical inoculating chain, fermentation process and equipment.

These two strains were evaluated for the production of vanillin by reducing to practice the process taught in US Patent 6,133.003 to Rabenhorst et al. asmetines referred to as the Symrise process in the following. At the end of the present memo, some results are given while using strain Zyl 926 (Asp Lii-2) according to the process taught in US Application № 10/550,945 to Heald et al, referred to as the Saltais process.

## 2. Material and Methods

#### 2.4 Cultures work out

The process applied in this study was according to the teaching of US Patent 6.133,003 to Rabenhorst et al "Process for the preparation of vanillin and microorganism suitable therefor".

#### 2.1.1 Inoculating chain

- According to working Example 1 of US 6,133,003

Preparation of the preliminary culture

A 500 ml conical flask with a side baffle was filled with 100 ml of medium, comprising 1 g of mait extract. 0,4 g of glucose and 0,4 g of yeast extract and made up to 100 ml with water and then sterilized for 20 min at 121°C. After cooling, the flask was inoculated with 200 pt of a frozen glycerol culture of Amycolatopsis sp. DSM 9992 or Zyt 926 (Asp L-ii 2). The cultures were incubated on a rotary shaking machine at 45 °C and 100 rpm. After 24 h, these cultures were used to inoculate the production medium.

- Preliminary culture medium preparation



| Product                    | Concentration |  |
|----------------------------|---------------|--|
| Glucose                    | 4 g/Kg        |  |
| Malt extract               | 10 g/Kg       |  |
| Yeast extract              | 4 g/Kg        |  |
| H <sub>2</sub> O q.s.p     |               |  |
| pH not adjusted            |               |  |
| Sterilization 20' at 121°C |               |  |

## · Preculture Inoculation

\*inoculation of 100 mL with 9,2 mL from the content of a cryotube (frozen glycerol culture)

"Incubation (shaking machine): 45°C - 100 rpm - 24h

## 2.1.2 Cultures in fermenters (maximum working volume = 15L)

Stage production according to working Example 3 of US 6,133,003

Vanillin Production in a 10 L fermenter.

5.1 of culture medium (4 g/l of glucose, 10 g/l of malt extract and 6 g/l of yeast extract) were sterilized in a fermeter and, after cooling, were inoculated with 100 ml of a seed culture of DSM 9992 or Zyl 926 (Asp L-ii 2) according to Example 1.

According to working Example 3, the culture conditions were:  $37^{\circ}$ C, 500 rpm, 51 of air/min, 12,5 h after inoculation, 1.634 kg of an approximately 3.7% strength ferulic acid solution (60.2 g of ferulic acid) were added. After 17.5 h, a further 4.397 kg of an approximately 3.7% strength ferulic acid solution (164.72 g of ferulic acid) were pumped in over a period of

While trying to reproduce these very same parameters, I had to adapt some of the above data to make it work and arrive at the following parameters:

LESAPERE



## - Culture parameters

|                   | GROWTH                                                                    | BIOCONVERSION                                                        | BIOCONVERSION |  |  |
|-------------------|---------------------------------------------------------------------------|----------------------------------------------------------------------|---------------|--|--|
| Inoculation       | 100 mL                                                                    | 100 mL                                                               |               |  |  |
| Temperature       | 37°C                                                                      | 37°C 37°C                                                            |               |  |  |
| Medium quantity   | 5 Kg                                                                      | 1.6 Kg 4.4 k                                                         |               |  |  |
| Culture condition | Batch                                                                     | Batch (isolated addition Fed-batch during of ferulic acid) (7.2 mL/m |               |  |  |
| Duration          | 12.5h                                                                     | 12.5h 5h 14.5i                                                       |               |  |  |
| pO₂ Regulation    | 50% (not mentioned in the example 3 of the patent)                        |                                                                      |               |  |  |
| Aeration          | 5 L/min                                                                   | 5 L/min 7 L/min 7 L/min                                              |               |  |  |
| Stirring          | 400 à 750 rpm                                                             |                                                                      |               |  |  |
| pH regulation     | $8\pm0.5$ (not mentioned in the example 3 of the patent but in the spec.) |                                                                      |               |  |  |
| Surpression       | No                                                                        |                                                                      |               |  |  |
| End of culture    | -                                                                         | - 32 h                                                               |               |  |  |

## Fermenters preparation

|                                 | ,             |  |
|---------------------------------|---------------|--|
| Product                         | Concentration |  |
| Glucose                         | 4 g/Kg        |  |
| Mait extract                    | 10 g/Kg       |  |
| Yeast extract                   | 6 g/Kg        |  |
| Propylene glycol<br>P2000       | 0.35 g/Kg     |  |
| H₂O                             | q.s.p         |  |
| 5 Kg of medium in the fermenter |               |  |

5 Kg of medium in the fermenter in situ sterilization



- Ferulic acid solution

| Product             | Concentration |
|---------------------|---------------|
| Ferulic acid        | 3.75%         |
| pH adjustment : 7.4 |               |

rentum annu hanta te dinembra fersila naid in umban yang namand nuh in umba. I kant ta maid

Detailed sequence of addition of the ferulic acid solution

- At T12,5 h of culture, 1,6 kg of the ferulic solution were added in the fermenter
- For the semi-continuous step which starts at 17,5 h of culture, 1,3 kg of solution were prepared to feed in the fermenter. If was necessary to reload the solution every 2,5 h because a degradation of fertilic and is observed when the solution is prepared more than 3 h in advance. The flow rate is 7,2 m/min during 10 h, allowing to add 4,4 kg of fertilic acid solution.
- After 32 h, the fermentation was terminated.

## 2.2 Culture follow up

## 2.2.1 Culture parameters

For each culture, parameters were monitored as follows:

- pO2 (dissolved oxygen in the medium) which allows to visualize the bacterial growth by the more or less important oxygen demand. This parameter is controlled by regulating aeration and/or agitation of the medium to a certain level in order to secure optimal growth conditions
- pH is regulated to favor bacterial growth and the production of metabolites of interest
- Solutions used for pH regulation [sulfuric acid (H<sub>2</sub>SO<sub>4</sub>) and/or sodium hydroxyde (NaOH)]

## 2.2 2 Bacteriological follow up

A bacterial follow up was realized in order to be ascertain the absence of any contamination.

## 2.2.3 HPLC Analyses

For each culture, several samples were withdrawn at mentioned times and frozen until analyzed:

<sup>\*</sup>preliminary tests to dissolve ferulic acid in water were carned out in valin, I had to acid ecidium hydroxide otherwise I could not get any homogeneous solution.



TO TA

TR

T12.5 after spot addition of ferulic acid

T17,5 after the starting of the semi-continuous culture

T20 before reloading of the ferulic acid solution

T22,5 before reloading of the ferulic acid solution

T25 before reloading of the ferulic acid solution

T27,5 before termination of the semi-continuous culture

Production of vanillin, vamilic acid, vanilic alcohol, gualacol and of the ferulic acid consumption were assayed by HPLC (high pressure liquid chromatography) and UV detection, consumption of glucose was tested through HPLC and retractometry

Two dilutions were made for each sample in order to be sure to be in the detection range of the metabolites.

| 3. | Results                                 |  |
|----|-----------------------------------------|--|
|    | *************************************** |  |

## 3.1 Culture follow up

## 3.1.1 DSM 9992

As illustrated by Figure 1, growth phase went fine, as shown by the pO2 which decreased from 100 % to 57 % between T0 and T12,5 h. This pO2 value had slightly risen during the spot addition of feruic acid (to 61 %) then decreased again regularly to reach 36 % at the termination of the culture. The 50 % instruction value for pO2 could not be maintained, though agitation was set at a high limit but pO2 did not drop under 30 %, thus oxygen supply has been sufficient throughout the culture duration.


Regarding pH regulation, the addition of acid and/or base solution was not necessary, let alone at the starting of the culture to restore the pH at  $8.0\pm0.5$ .

At T24 h, I observed a high toam production necessitating the addition of antifoam propylene glycol P2000.

At the termination of the culture, the culture medium was of a yellowlorange color, greasy —probably due to the presence of artifloam) and I noticed very high levels of clove ador/spice, showing a priori that bioconversion practically did not occur.











Figure 1: Parameter acquisition for DSM 9992 strain cultivated according to Symrise process (US 6,133,033) (Cumul = Total)


## 3.1.2 Zvl 926 (Asp L-ii 2)

As shown in Figure 2, this experiment took place with noticeable differences compared to the previous one with strain DSM 9992:

- pO2 decreased regularly during growth phase then went over 80 % after the spot addition of ferulic acid
- pH regulation by base solution addition only after addition of ferutic acid during the bioconversion (total amount added: 120g of NaOH 10%).










Figure 2 : Parameter acquisition for the Zyl 926 (Asp L-ii 2) strain cultivated according to the Symrise process (US 6,133,033)

## 3,1,3 Cultures comparison

Table 1 illustrates the principal data:






Table 1: Growth comparison of DSM 9992 and Zvl 926 (Asp L-ii 2) strains

|                         | pO2                                                                                                                                                                                        | рН                                                                                                         | pH regulation<br>with sodium<br>hydroxide 10%                                                                 | pH regulation<br>with sulfuric acid<br>10% | Medium<br>composition at<br>the end of the<br>culture                                  |
|-------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|--------------------------------------------|----------------------------------------------------------------------------------------|
| DSM 9992                | Growth phase (between 0 - 12 h): decrease from 100% to 57%  Bioconversion phase (between 12 - 32 h): increase to 61% after ferulic acid addition, then decrease until the end (35% at 32h) | Low initial pH $(5.4)$ $\rightarrow$ addition of 15g of sodium hydroxide 10% to adjust pH at $8.0 \pm 0.5$ | Bioconversion<br>phase (between<br>12 - 32 h) : NO                                                            | NO                                         | Oily medium<br>(due to<br>antifoam),<br>yellow/orange,<br>strong smell like<br>a clove |
| Zyl 926<br>(Asp L-ii 2) | Growth phase (between 0 - 12 h): continuous decrease from 100% to 69%  Bioconversion phase (between 12 - 32 h): increase to 85% after ferulic acid addition, then static until the end     | No need to adjust initial pH                                                                               | Continuous<br>addition of<br>sodium<br>hydroxide 10%<br>during<br>bioconversion<br>phase<br>(quantity : 120g) | NO                                         | Orange<br>medium, smell<br>like butter                                                 |

## 3.2 Analytical results

Table 2 allows to compare the results after the termination of the culture for the two strains :

Table 2: End product concentrations for Zyl 926 (Asp L-ii 2) and DSM 9992 strains when cultivated according to Example 3 of US 6,133,003

| Strain     | Vanillic<br>alcohol<br>(g/Kg) | Vanillic acid<br>(g/Kg) | Vanillin<br>(g/Kg) | Guaiacol<br>(g/Kg) |
|------------|-------------------------------|-------------------------|--------------------|--------------------|
| Asp L-ii 2 | 0                             | 0                       | 0                  | 0                  |
| DSM 9992   | <0.1                          | 0.3                     | <0.1               | 0                  |



Chromatograms below clearly indicate that there was an accumulation of feruitic acid over time, whatever strain was used.

No metabolite of interest was detected for strain Zyl 926 (Asp L-ii 2) as illustrated in Figure 3 below. and as shown in Figure 4 below, almost none in the case of strain DSM 9992 : under 0.1 g/kg (0.1 g/l) regarding vanillin, which value should be interpreted as follows : practically no vanillin was produced contrary to the statement in Example 3 of US 6,133,003.

Moreover at T17.5 h-T20 h-T22.5 h-T25 h and T27.5 h analysis of intermediary samples was conducted and showed no evidence of the presence of metabolites which could have been produced and then eventually degraded during the cultivation period.




Figure 3: Chromatogram/Symrise process (US 6.133033) with Zvl 926 strain (Asp L-ii 2) sample corresponding to the end of the culture (T32 h) vaniline/vanilin -- guaiacol/guaiacol -- acide férul-que/ferul-c acid)




Figure 4 : Chromatogram/Symrise process (US 6,133,033) with DSM 9992 strain – sample corresponding to the end of the culture (T32 h)

[French/English translation: alcool vanifique/vanific alcohol – acide vanifique/vanifilic acid – vanifiline/vanifin – guaracol/guaracol – acide férulique/ferulic acid

Table 3 illustrates the accumulation profile of ferulic acid during cultures.

Table 3: Accumulation profile of ferulic acid in the culture medium

| Sample                                           | Ferulic acid<br>Theoretical<br>quantity added (g) | Ferulic acid<br>Theoretical<br>quantity<br>cumulated* (g) | Ferulic acid<br>Theoretical<br>concentration*<br>(g/Kg) | Ferulic acid<br>HPLC results<br>(g/Kg) |
|--------------------------------------------------|---------------------------------------------------|-----------------------------------------------------------|---------------------------------------------------------|----------------------------------------|
| T12.5 after addition of<br>ferulic acid solution | 60.2 g                                            | /a                                                        | 9.1 g/Kg                                                | 8.9 (± 0.3)**                          |
| T20 before fed-batch<br>reloading                | 48.75 g                                           | 108.9 g                                                   | 13.8 g/kg                                               | 15.0 (± 2.4)**                         |
| T22.5 before fed-batch<br>reloading              | 48.75 g                                           | 157.7 g                                                   | 17.1 g/Kg                                               | 17.3 (± 1.8)**                         |
| T25 before fed-batch<br>reloading                | 48.75 q                                           | 206.5 q                                                   | 19 7 g/Kg                                               | 18.5 (± 1.6)**                         |
| 127.5 at the end of fed-<br>batch culture        | 48.75 g                                           | 255.3 g                                                   | 21.6 g/Kg                                               | 21.1 (± 0.8)**                         |
| T32                                              |                                                   |                                                           |                                                         | 19.8 (± 1.1)**                         |

<sup>&</sup>quot;Two hypothesis

<sup>-</sup> no addition of sodium hydroxide to regulate pH

<sup>-</sup> no consumption and/or loss of ferulic acid between fed-batch reloading

<sup>\*\*</sup> Result average (g/Kg)



As illustrated by theoretical calculations reported in Table 3 above, bioconversion phase was conducted as taught in US patent 6.133,003 but feruito acid accumulated and its consumption and degradation after introduction in the culture medium are negligible.

#### 4. Conclusion

Though conducted with identical culture conditions, the above analytical results allow to conclude that strain DSM 9992, claimed in US 6,133,003 and 2yl 926 (Asp. L-ii 2), claimed in US Application N° 10/SS0,945, have shown a different behavior when grown in fermenters. Besides, in both cases growth phase went well.

Moreover, HPLC analysis were performed in order to quantify metabolities of interest, namely vanillin, guaracot, vanillin alcohol and vanillic acid expected to be produced in such conditions.

These trials results allow us to conclude that no ferulic acid reacted and the latter accumulated in the culture medium during the course of the culture.

These were unexpected results since the parameters applied in the present study were those taught in working. Examples of the document US. 6,133,003 to Rabenhorst et al., "Process for the preparation of vanillin and microorganisms suitable therefor".

Moreover we, present engineers who signed the appended affidavite, made all reasonable attempts to use the complete teaching of this specification and not only the Example wording. See for example to name but a lew "adjustments", the preparation of feruitic acid and also pH and pO<sub>2</sub> regulation and antifoam necessary usage.

As a conclusion, and based on the forgoing facts, this prior art, namely US 6,133.903 is of a nonoperative kind, the reason being that a person skill in the art by repeating/reproducing its teaching not only does not get the expected results but rather gets essentially no production of the desired product, namely vanillin. To arrive at this production, one should "invent" another process.

Accordingly, US 6,133,003 is not at all pertinent, and by far, the production of vanillin as taught and claimed by the inventors of US Application N° 10/650,945 is not obviously derivable therefrom.

By way of comparison below is shown in Figure 5, a chromatogram obtained by working out the process faught in US Application N° 10/550,945 to Heald et al., namely the "Safisis process", using Zyl 926 (Asp L-ii 2): fendic acid was spent and the vanillin peak is predominant (14 g/l).






Figure 5 : Chromatogram/Safisis process (US App. N° 10/550,945) with Zyl 926 strain (Asp L-ii 2) - sample withdrawn at the end of the culture (French/English translation, abod vinificial/white labell - acide valificial/white in - guitaculture in - guitaculture - acide (fruit/qualficulture) - acide (fruit/